

Measure the fluorescence excitation and emission spectrum

1. Turn on the instrument
 - a. Turn on the lamp power
 - b. Turn on the computer, monochromator, detector.
 - c. Open the software Felix
2. Excitation measurement procedure
 - a. In the menu 'Acquisition', select 'New Acquisition/excitation scan'
 - b. Choose the appropriate excitation wavelength range and the emission wavelength. (The emission and excitation wavelength cannot overlap)
 - c. Click the button 'More', select 'Real time spectrum correction' and select both 'Excitation correction' and 'emission correction'
 - d. Click 'Acquire', then 'Start'
3. Emission measurement procedure
 - a. In the menu 'Acquisition', select 'New Acquisition/Emission scan'
 - b. Choose the appropriate emission wavelength range and the excitation wavelength. (The emission and excitation wavelength cannot overlap)
 - c. Click the button 'More', select 'Real time spectrum correction' and select both 'Excitation correction' and 'emission correction'
 - d. Click 'Acquire', then 'Start'
4. Measure Fluorescein excitation and emission spectrum
 - a. Measure the buffer excitation spectrum (excitation: 400nm-510nm, emission: 530nm) and emission spectrum (excitation: 460nm, emission: 480nm-600nm)
 - b. Dilute the fluorescein stock solution by a factor 100 x 50 (16nM)
 - c. Measure the excitation and emission spectra using the same settings as in a.
 - d. Subtract the measured fluorescein spectrum from the background spectrum, normalize; you now have the corrected spectrum.
 - e. Now take the emission spectrum at a different excitation wavelength (440nm). Compare the emission spectra at different excitation wavelengths.
 - f. Dilute the fluorescein stock solution by a factor of 100x100 (8nM), measure the emission spectrum. Compare the fluorescence spectra and show that the concentration and fluorescence signal are proportional to each other.
5. Measure Fluorescence anisotropy
 - a. Put the polarizer in the excitation and emission path
 - b. In the 'Acquisition', select 'New Acquisition/Emission Scan with Polarizer'. Choose appropriate excitation and emission wavelengths for fluorescein.
 - c. Calibrate the G-factor: Use horizontal excitation light, measure vertical and horizontal emission fluorescence. Take the ratio between vertical and horizontal fluorescence.

- d. Measure the anisotropy of the fluorescein sample. Use vertical excitation light, measure the vertical and horizontal emission fluorescence. Calculate the anisotropy.
- e. Add 4420 antibody and measure the anisotropy.